

ASSIGNMENT NO: 4 : EXPERIMENT ON Ph & CONDUCTIVITY MEASUREMENT

❖ SAMPLE: LEMON JUICE

Level-1 Calibration

Calibration Diagram:

Control Panel:

Calibration at 25°C

Initial pH Value:

pH Value:

pH	OUTPUT mV
0	414
4	177
7	0
14	-414

Plot **Reload** **Level-2**

Level-2 Measurement

Measurement Diagram:

Control Panel:

Temperature at 25°C

Sample:

Enter Output Voltage:
(in mV)

Submit

Sample	pH	OUTPUT mV
Lemon Juice	2.2	283.97

Reload **Level-3**

Level-3 Temperature

[Temperature Diagram:](#)

High Impedance Buffer Amplifier Display

[Control Panel:](#)


Sample: Lemon Juice
At 25°C Output mV = 283.97



[Get Temperature =>](#) 13°C

Temperature °C	pH	OUTPUT mV
13	2.2	272.54

[Plot](#) [Reload](#)

[Formula](#)





❖ SAMPLE : COFFEE

Level-2 Measurement

Measurement Diagram:

High
Impedance
Buffer Amplifier Display

Control Panel:

Temperature at 25°C

Sample:

Enter Output Voltage:
(in mV)

Submit

Sample pH OUTPUT mV

Coffee 5 118.32

Reload Level-3

<--Level-1

Sensor Analysis Laboratory

pH Calibration

Level-3 Temperature

Temperature Diagram:

High
Impedance
Buffer Amplifier Display

Control Panel:

Sample:

At 25°C Output mV = 118.32

Get Temperature → 36°C


Temperature °C pH OUTPUT mV

36 5 122.68

Plot Reload

<--Level-2

Formula



❖ SAMPLE : AMMONIA

Level-2 Measurement

[<-Level-1](#)

Measurement Diagram:

High
Impedance
Buffer Amplifier Display

Control Panel:

Temperature at 25°C

Sample:

Enter Output Voltage:
(in mV)

Submit

Sample	pH	OUTPUT mV
Coffee	5	118.32
Ammonia	11	-236.64

Buttons: Reload, Level-3, Formula

Level-3 Temperature

[<-Level-2](#)

Temperature Diagram:

High
Impedance Amplifier Display
Buffer

Control Panel:


Sample:

At 25°C Output mV = -236.64

Get Temperature ==> 24°C

Temperature °C	pH	OUTPUT mV
24	11	-235.84

Buttons: Plot, Reload, Formula



❖ **SAMPLE: NaOH**

Sensor Analysis Laboratory

Measurement of Conductivity

Level-1 Measurement

Control Panel:

- Sample: NaOH
- Concentration: 30%
- Cell Constant: 1.0
- Specific Conductance at 25°C: 6.17 (Siemens)

Console window:

To Have a better Comparison Selected Sample will remain same for next two levels

Selected Substance is NaOH
Selected Concentration is 30
Your Answer : 6.17

Buttons: Plot, Submit, Next Set of Values, Level-2, Reload, Formula

Right Panel:

- PPM Concentration : 300000
- Equivalent weight : 40.01
- Density at 25 °C : 1.331

Level-2 Temperature

Control Panel:

- Sample: NaOH
- Specific Conductance at 25°C = 6.17
- Concentration: 30%

Get Temperature => 38°C

Conductivity at Temperature of Interest

7.78 (Siemens)

Buttons: Plot, Submit, Next Set of Values, Level-3, Reload, Formula

Right Panel:

- Temperature Range from: 0°C to 50°C
- Temperature Co-efficient : 0.0201

Level-3 Contamination

←Level-2

Control Panel:

Sample: NaOH

Concentration: 30%

Cellconstant: 1.0

Contamination: 0.2mm

Temperature Range from: 0°C to 30°C

Temperature Co-efficient : 0.0201

Formulas

Console window:

Assumption: Half of this value is deposited on each electrode.

Default values for L, A and B:
L=1cm , A=10cm , B=0.1cm

Selected Contamination is 0.2

Modified Cell Constant : 0.82

Specific Conductance at 25°C : 6.17 (Siemens)

Modified specific conductance value: 5.06 (Siemens)

Reload

❖ SAMPLE: H₂SO₄

Level-1 Measurement

Control Panel:

Sample: H₂SO₄

Concentration: 50%

Cell Constant: 1.0

Specific Conductance at 25°C :

8.06 (Siemens)

Plot Submit

Next Set of Values Level-2

Reload

Formula

PPM Concentration : 500000

Equivalent weight : 87.14

Density at 25 °C : 1.399

Console window:

To Have a better Comparison Selected Sample will remain same for next two levels

Selected Substance is NaOH
Selected Concentration is 30
Your Answer : 6.17
Selected Substance is H₂SO₄
Selected Concentration is 50.
Your Answer : 8.06

Level-2 Temperature

Control Panel:

Sample: H₂SO₄

Specific Conductance at 25°C = 8.06

Concentration: 50%

Get Temperature => 33°C

Conductivity at Temperature of Interest

9.30 (Siemens)

Plot Submit

Next Set of Values Level-3

Reload

Formula


Temperature Range from: 0°C to 50°C

Temperature Co-efficient : 0.0193

Console window:

Selected Concentration is 30
Your Answer ; 7.78
Selected Concentration is 50
Your Answer : 9.3

Level-3 Contamination



Control Panel:

Sample: H2SO4

Concentration: 50%

Cellconstant: 1.0

Contamination: 0.2mm

Modified Cell Constant: 0.82

Specific Conductance at 25°C: 8.06 (Siemens)

Modified specific conductance value: 7.8 (Siemens)

Console window:

Assumption: Half of this value is deposited on each electrode.

Default values for L, A and B:
L=1cm , A=10cm , B=0.1cm

Selected Contamination is 0.2

Temperature Range from: 0°C to 50°C

Temperature Co-efficient : 0.0193

Formula

Reload

❖ SAMPLE: HCL

Control Panel:

Sample: HCL

Concentration: 10%

Cell Constant: 1.0

Specific Conductance at 25°C:

1.79 (Siemens)

Plot Submit


Next Set of Values Level-2

Reload

Formula

PPM Concentration : 100000
Equivalent weight : 36.46
Density at 25 °C : 1.048

Console window:
To Have a better Comparison Selected Sample will remain same for next two levels
Selected Substance is NaOH
Selected Concentration is 30
Your Answer : 6.17
Selected Substance is H₂SO₄
Selected Concentration is 50
Your Answer : 8.06
Selected Substance is HCl
Selected Concentration is 10
Your Answer : 1.79



Level-2 Temperature

Control Panel:

Sample: HCL

Specific Conductance at 25°C = 1.79

Concentration: 10%

Get Temperature ==> 31°C

Conductivity at Temperature of Interest

1.96 (Siemens)

Plot Submit


Next Set of Values Level-3

Reload


Formula

Temperature Range from: 0°C to 50°C
Temperature Co-efficient : 0.0156

Console window:
Selected Concentration is 30
Your Answer : 7.78
Selected Concentration is 50
Your Answer : 9.3
Selected Concentration is 10
Your Answer : 1.96



Level-3 Contamination



Control Panel:

<-Level-2

Temperature Range from: 0°C to 50°C

Temperature Co-efficient : 0.0156

Console window:

Assumption: Half of this value is deposited on each electrode.

Default values for L, A and B:
L=1cm , A=10cm , B=0.1cm

Selected Contamination is 0.2

Sample: HCL

Concentration: 10%

Cellconstant: 1.0

Contamination: 0.2mm

Modified Cell Constant : 0.82

Specific Conductance at 25°C : 1.79 (Siemens)

Modified specific conductance value: 1.68 (Siemens)

Reload

Formula

142103005

S.Y. Comp

Sensors and AutomationAssignment No: 4

Q1. What is the effect of temperature measurement of pH?

Answer When the temperature of a solution rises.

The modulus vibrations in the Solution rise resulting in the ionization and formation of H^+ ions. More H^+ ions lead to more acidic behaviour. During the temperature changes, the pH value of the solution changes. Thus pH decreases on increasing the temperature.

Q2. Describe the construction and principle of working of a conductivity meter.

Answer Construction :

The meter is equipped with a probe, usually handheld for field or on-site measurement. After the probe is placed in the liquid to be measured the meter applies voltage between two electrodes inside the probe. Electrical resistance from the solution causes a drop in voltage, which is read by the meter. The meter converts this reading to milli or microsiemens per centimeter. The value indicates the total dissolved solids. Total dissolved solids is the amount of solids that can pass through a glass fiber filter.

Principle of working:

Two plates are placed in the sample, a potential is applied across the plates and the current that passes through the solution is measured.

3. Describe the calibration procedure for the pH meter.

Answer 1. The calibration must be performed for buffers with pH 4.0, 7.0 and 9.2

At least a 2-point calibration must be performed at room temperature using buffers that meet the expected pH value of the sample.

If a one point calibration is performed, measurement errors are more for the sample that is being measured.

If the temperature probe is not used at the time of calibration, default value considered is 25°C .

4. Discuss any three applications where the conductivity meter is used.

Answer 1. Hydroponics:

The higher the reading on conductivity meter, the more nutrients that are available for your plants. But it's easy to over-fertilize the indoor garden and end up with nutrient burn so the best EC levels are in the moderate range during the growing period about 1.2 to 1.6 and no longer or higher than 1.8 during flowering.

2. Purity of water:

Conductivity measures water's ability to conduct electricity due to presence of / or absence of certain ions. While pure water conducts electricity poorly, water that has certain chemicals or elements in it and at varying amounts is a better conductor of electricity.

142103005, SY Comp

3. Aquaculture:

For water resources, the conductivity may indicate whether or not the water is too saline to be drinkable or usable for irrigation or industrial use.

5. What are the advantages and disadvantages of a conductivity meter?Answer Advantages:

Conductivity offers a fast, reliable, nondestructive, inexpensive and durable means of measuring the ionic contents of a sample. Reliability and repeatability are excellent.

Disadvantages:

1. It cannot distinguish between different types of ions.
2. It gives instead a reading proportional to the combined effect of all ions present.
3. It must be applied with some pre-knowledge of the solution composition or used in relatively pure solutions to be successful.

142103005

- Calculations:

- Conductivity at temp of interest:

$$T = 38^\circ\text{C} = 38 + 273 = 311\text{K}$$

$$CR = 6.17, \quad t_C = 0.02011$$

$$C_T = CR(1 + C(T-25)) \text{ Siemens}$$

$$C_T = 6.17 (1 + 0.0201 (38 - 25))$$

$$\therefore \underline{C_T = 7.78 \text{ Siemens}}$$

- Lemon Juice

$$T = 25^\circ\text{C} = 298\text{K}$$

$$E_0 = 0, \quad F = 96485, \quad R = 8.314 \text{ J mol/K}$$

$$pH = 2.2; \quad pH_e = 2.2 - T = -4.8$$

$$E = E_0 - 2.303 \left(\left(\frac{R \times T}{F} \right) \times pH_e \right) \times 108$$

$$\underline{E = 283.97 \text{ mV.}}$$

- Coffee

$$T = 25^\circ\text{C} = 298\text{K} \quad E_0 = 0, \quad R = 8.314 \text{ J mol/K}$$

$$F = 96485 \text{ C/mole}$$

$$pH = 5, \quad pH_e = 5 - 7 = 2$$

$$E = E_0 - 2.303 \left(\left(\frac{R \times T}{F} \right) \times pH_e \right) \times 1000$$

$$\underline{E = 118.25 \text{ mV}}$$

- Ammonia:

$$T = 25^\circ\text{C} = 298\text{K}; \quad E_0 = 0, \quad R = 8.314 \text{ J mol/K}$$

$$F = 96485 \text{ C/mole}$$

$$pH = 11, \quad pH_e = 11 - 7 = 4$$

$$E = E_0 - 2.303 \left(\left(\frac{R \times T}{F} \right) \times pH_e \right) \times 1000$$

$$\underline{E = -236.63 \text{ mV.}}$$

142103005

• NaOH

concentration = 30%

$$\text{Specific conductance} = 1000 C A_0 (1 - a\sqrt{C} + bC)$$

$$\text{Here, } A_0 = 246.5,$$

$$a = 0.47$$

$$b = 0.3$$

$$T = 25^\circ\text{C} = 298\text{K}$$

$$\text{Normality} = \frac{3 \times 10^5 \times 1.33}{10^3 \times 40.01}$$

$$= 9.97 \times 10^6$$

$$\therefore \text{Specific Conductance} = \frac{1000 \times 9.8 \times 246.5}{(1 - 0.41\sqrt{9.97} + 5.3 \times 9.97)}$$

$$= 24156 \times 10^6 (1 - 4.483 + 0.34997)$$

$$CR = 6.17 \text{ Siemens.}$$

• H₂SO₄

Concentration = 50% , $A_0 = 151.4$ $a = 1.24$ $b = 1.14$

$$C(\text{Normality}) = \frac{\text{ppm l}}{1000 \times \text{eq. wt}} = \frac{5 \times 1.399 \times 10^5}{1000 \times 97.14}$$

$$C = 8.027$$

$$\therefore \text{Specific conductance} = 1000 A_0 (1 - a\sqrt{C} + bC) \text{ siemens}$$

$$= 1000 \times 8.027 \times 151.4 (1 - 1.24\sqrt{8027} + 1.14 \times 8027)$$

$$= 8.069 \text{ siemens.}$$

• Conductivity of at temperature of interest

$$T = 33^\circ\text{C} , C_p = 8.06 , t_c = 0.0193$$

$$C_T = CR [1 + t_c (T - 25)] \text{ siemens}$$

$$C_T = 8.069 [1 + 0.0193 (33 - 25)]$$

$$C_T = 9.3 \text{ siemens.}$$

142103005

• HCl

$$\text{concentration} = 10\%$$

$$A_0 = 420, a = 0.37, b = 0.38$$

$$C(\text{Normality}) = \frac{10^5 \times 1.048}{36.46 \times 10^3}$$

$$= 2.87$$

$$\therefore \text{specific conductance} = 1000 A_0 C (1 - aJ_C + b\chi)$$

$$= 1.79 \text{ siemens}$$

conductivity at temperature of interest

$$T = 31^\circ\text{C},$$

$$C_R = 1.4, T_C = 0.0156$$

$$C_T = C_R [1 + t_C (T - 25)] \text{ siemens}$$

$$C_T = 1.791 + 0.056 [31 - 25]$$

$$C_R = 1.957 \text{ siemens.}$$